

APPROXIMATE WEIGHT IN POUNDS

	COST CRUNCHER	TRUELITE	STARLITE	PRESTIGE
.50MM WIDE BELT	61420	66400	46400	56400
Small	7.4	6.8	6.2	5.6
Medium	8.3	7.7	7.0	6.3
Large	9.6	8.8	8.0	7.2
XL	10.9	10.0	9.1	8.2
XXL	12.4	11.4	10.3	9.3
.35MM WIDE BELT	31420	36400	86400	76400
Small	5.5	5.1	4.7	4.2
Medium	6.2	5.7	5.2	4.7
Large	7.1	6.5	6.0	5.4
XL	8.0	7.4	6.7	6.1
XXL	9.1	8.4	7.6	6.9
.50MM DELUXE	60947	67400	47400	57400
Small	7.3	6.7	6.1	5.5
Medium	8.2	7.5	6.8	6.1
Large	9.5	8.7	7.9	7.1
XL	10.7	9.8	8.9	8.0
XXL	12.2	11.2	10.1	9.1
.35MM DELUXE	30947	37400	87400	77400
Small	5.3	4.9	4.5	4.1
Medium	6.0	5.5	5.0	4.6
Large	6.9	6.4	5.8	5.2
XL	7.8	7.2	6.6	5.9
XXL	8.9	8.2	7.5	6.7
.50MM/.25MM STD FEMALE VEST WITH VELCRO	V60605	V66605	V46605	V56605
XS	4.7	4.3	3.9	3.5
Small	5.1	4.7	4.3	3.8
Medium	5.5	5.1	4.6	4.1
Large	5.9	5.4	4.9	4.4
XL	6.5	5.9	5.4	4.8
2XL	6.8	6.3	5.7	5.1
3XL	7.2	6.6	6.0	5.4
.40MM/.20MM FEMALE VEST WITH VELCRO	V30605	V36605	V86605	V76605
XS	3.9	3.6	3.2	2.9
Small	4.2	3.9	3.5	3.2
Medium	4.5	4.2	3.8	3.4
Large	4.8	4.4	4.0	3.6
XL	5.3	4.9	4.4	4.0
2XL	5.6	5.1	4.7	4.2
3XL	5.9	5.4	4.9	4.4
.50MM/.25MM FEMALE VEST WITH BUCKLES	V60321	V64321	V44321	V54321
XS	3.9	3.5	3.2	2.9
Small	4.2	3.9	3.5	3.2
Medium	4.6	4.2	3.8	3.5
Large	4.8	4.4	4.0	3.6
XL	5.8	5.3	4.8	4.3
2XL	6.2	5.7	5.1	4.6
3XL	6.6	6.1	5.5	5.0
.40MM/.20MM FEMALE VEST WITH BUCKLES	V30321	V34321	V84321	V74321
XS	3.2	2.9	2.6	2.4
Small	3.5	3.2	2.9	2.6
Medium	3.8	3.5	3.2	2.9
Large	4.0	3.6	3.3	3.0
XL	4.7	4.3	4.0	3.6
2XL	5.1	4.7	4.2	3.8
3XL	5.5	5.0	4.5	4.1
.50MM/.25MM REVERSE FEMALE VEST	V60605R	V66605R	V46605R	V56605R
XS	6.3	5.8	5.2	4.7
Small	6.9	6.3	5.7	5.1
Medium	7.4	6.7	6.1	5.4
Large	7.8	7.1	6.4	5.8
XL	8.3	7.6	6.8	6.1
2XL	8.8	8.1	7.3	6.5
3XL	9.3	8.5	7.7	6.9
.40MM/.20MM REVERSE FEMALE VEST	V30605R	V36605R	V86605R	V76605R
XS	5.2	4.7	4.3	3.8
Small	5.6	5.2	4.7	4.2
Medium	6.0	5.5	5.0	4.5
Large	6.4	5.8	5.3	4.7
XL	6.7	6.2	5.6	5.0
2XL	7.2	6.6	6.0	5.3
3XL	7.6	6.9	6.3	5.6

APPROXIMATE WEIGHT IN POUNDS

	COST CRUNCHER	TRUELITE	STARLITE	PRESTIGE
.50MM/.25MMSTD SKIRT UNISEX	S60862	S64862	S44862	S54862
XS	5.5	5.1	4.6	4.2
Small	5.9	5.4	4.9	4.5
Medium	6.5	6.0	5.4	4.9
Medium 2	6.8	6.3	5.7	5.2
Large	7.1	6.5	5.9	5.4
Large 2	7.7	7.1	6.4	5.8
XL	8.5	7.8	7.0	6.4
2XL	8.9	8.2	7.4	6.7
3XL	9.7	8.9	8.1	7.3
.40MM/.20MM SKIRT UNISEX	S30862	S34862	S84862	S74862
XS	4.5	4.2	3.8	3.5
Small	4.9	4.5	4.1	3.7
Medium	5.4	4.9	4.5	4.1
Medium 2	5.6	5.2	4.7	4.3
Large	5.9	5.4	4.9	4.4
Large 2	6.4	5.8	5.3	4.8
XL	7.0	6.4	5.8	5.3
2XL	7.3	6.7	6.1	5.5
3XL	8.0	7.3	6.7	6.0
.50MM/.25MM STD MALE VEST WITH VELCRO	V60650	V66650	V46605	V56605
XS	5.6	5.2	4.7	4.2
Small	6.1	5.6	5.0	4.5
Medium	6.9	6.3	5.7	5.2
Large	7.3	6.7	6.1	5.5
XL	7.9	7.2	6.6	5.9
2XL	8.1	7.4	6.7	6.0
3XL	8.6	7.9	7.1	6.4
.40MM/.20MM MALE VEST WITH VELCRO	V30605	V36605	V86605	V76605
XS	4.6	4.2	3.8	3.5
Small	5.0	4.6	4.1	3.7
Medium	5.7	5.2	4.7	4.3
Large	6.0	5.5	5.0	4.5
XL	6.5	5.9	5.4	4.8
2XL	6.6	6.0	5.5	4.9
3XL	7.0	6.5	5.9	5.3
.50MM/.25MM MALE VEST WITH BUCKLES	V60234	V61234	V41324	V51234
XS	4.6	4.2	3.8	3.4
Small	5.0	4.6	4.1	3.7
Medium	6.2	5.7	5.1	4.6
Large	6.6	6.0	5.5	4.9
XL	7.0	6.4	5.8	5.2
2XL	7.7	7.0	6.4	5.7
3XL	7.9	7.3	6.6	5.9
.40MM/.20MM MALE VEST WITH BUCKLES	V30234	V31234	V81234	V71234
XS	3.8	3.4	3.1	2.8
Small	4.1	3.7	3.4	3.1
Medium	5.1	4.7	4.2	3.8
Large	5.4	5.0	4.5	4.1
XL	5.8	5.3	4.8	4.3
2XL	6.3	5.8	5.2	4.7
3XL	6.5	5.9	5.4	4.9
.50MM/.25MM REVERSE MALE VEST	V60650R	V66650R	V46605R	V56605R
XS	7.6	7.0	6.3	5.6
Small	8.2	7.5	6.8	6.1
Medium	8.6	7.9	7.1	6.4
Large	9.3	8.5	7.7	6.9
XL	9.8	8.9	8.1	7.2
2XL	10.6	9.7	8.7	7.8
3XL	11.3	10.3	9.3	8.3
.40MM/.20MM REVERSE MALE VEST	V30650R	V36605R	V86605R	V76650R
XS	6.2	5.7	5.2	4.6
Small	6.7	6.1	5.5	5.0
Medium	7.0	6.4	5.8	5.2
Large	7.6	6.9	6.3	5.6
XL	8.0	7.3	6.6	5.9
2XL	8.6	7.9	7.1	6.4
3XL	9.2	8.4	7.6	6.8

APPROXIMATE WEIGHT IN POUNDS

	COST CRUNCHER	TRUELITE	STARLITE	PRESTIGE
.50MM/.25MM Wide Belt Wrap	60590	66590	46590	56590
Small	9.0	8.3	7.6	6.8
Medium	11.0	10.1	9.2	8.3
Large	12.3	11.3	10.3	9.3
XL	13.8	12.7	11.5	10.4
.40MM/.20MM Wide Belt Wrap	30590	36590	86590	76590
Small	7.5	6.9	6.3	5.7
Medium	9.1	8.3	7.6	6.9
Large	10.1	9.3	8.5	7.7
XL	11.4	10.5	9.5	8.6
.50MM/.25MM Comfort Wrap	60565	66565	46565	46565
Small	9.1	8.4	7.6	6.9
Medium	10.8	10.0	9.1	8.2
Large	12.5	11.5	10.4	9.4
XL	14.0	12.8	11.6	10.5
.40MM/.20MM Comfort Wrap	30565	36565	86565	76565
Small	7.6	7.0	6.4	5.8
Medium	9.0	8.2	7.5	6.8
Large	10.3	9.5	8.7	7.8
XL	11.5	10.6	9.6	8.7
.50MM/.25MM Stretch Back Reverse Wrap	60582	66582	46582	56582
Small	10.8	9.9	8.9	8.0
Medium	12.6	11.5	10.4	9.3
Large	14.3	13.1	11.8	10.5
XL	15.8	14.4	13.0	11.6
.40MM/.20MM Stretch Back Reverse Wrap	30582	36582	86582	76582
Small	8.8	8.1	7.3	6.5
Medium	10.2	9.3	8.5	7.6
Large	11.6	10.6	9.6	8.6
XL	12.8	11.7	10.6	9.5